




With neQis Sterilizer Safety Detection, you're in safe hands...

neQis Stops the Release of Non-Sterile Loads

In the working environments of today's sterile services departments, it is still possible that an operator can release a non-sterile load.

In some instances, non-sterile loads have only been highlighted after some items of a load have been used.

It can all happen by the operator placing the load within the sterilizer and then either being distracted or having to leave the department for a reason. On returning the operator opens the door and releases the load, not realising the sterilizer had not been run. It is not until the sets are to be used on a patient, that in some cases it is noticed that the tape & indicators have not changed colour.

Have you released an unsterile load?

Has a load been placed in an autoclave but not started?

'But it was hot when I removed it from the machine'

'I only got distracted for a moment'.....

Sound familiar?

The thought was then raised, if this one load has been caught, how many others haven't?

neQis have now developed a safe and effective system ensuring any load placed into a sterilizer has to be processed before the load can be removed.

Would you like to

Reduce the risk, install confidence, & save potential litigation. Now you can.

The system offers

- Door interlocks, ensuring the Sterilizer doors cannot be opened unless the machine has completed a cycle
- Fits to any manufacturers machine.
- Local alarm indicators
- Safety override Key- switch
- A major step towards patient safety.

To find out more call our team on

07809 702356 or **07738 111 355**

Email: info@neqis.co.uk

neQis Ltd

Independent Monitoring Systems